

Testing for the Differential Expression of Genes at the Probe Level of Affymetrix Microarray data

T. Khamiakova

Interuniversity Institute for Biostatistics and Statistical bioinformatics



Universiteit Hasselt

29 September 2010

Overview

Probe Level
Analysis

Khamiakova et
al.

Introduction

Affymetrix
Microarray
Probe Level
Data

Methods

LMM

Dataset

Results and
Discussion

Simulation
Study

1 Introduction

2 Affymetrix Microarray Probe Level Data

3 Methods

- Mixed model approach
- Dataset

4 Results and Discussion

5 Simulation Study

Advantage of Probe Level Analysis

Probe Level
Analysis

Khamiakova et
al.

Introduction

Affymetrix
Microarray
Probe Level
Data

Methods

LMM

Dataset

Results and
Discussion

Simulation
Study

Standard practice:

- summarization (with previous pre-processing)
- analysis of differential expression

Drawback of Summarization:

- Loss of information

Alternative:

- Differential expression on probe level

Probes:

- Probeset - multivariate measurement of gene expression
- Probe = random variable measuring expression of gene, has a specific mean μ_j

Data structure

- each array can add extra variability to GE measurements
- in more complex designs fixed effects also introduce variability

Illustration

Probe Level
Analysis

Khiamiakova et
al.

Introduction

Affymetrix
Microarray
Probe Level
Data


Methods

LMM

Dataset

Results and
Discussion

Simulation
Study



Adjustment for covariates

Probe Level
Analysis

Khamiakova et
al.

Introduction

Affymetrix
Microarray
Probe Level
Data

Methods

LMM

Dataset

Results and
Discussion

Simulation
Study

LMM-2:

$$\log(PM_{ijk}) = \mu_j + \alpha \cdot T_{ik} + b_i + \epsilon_{ijk}$$

- μ_j - probe-specific mean
- $b_i \sim N(0, \sigma_b^2)$ - array-to-array variability
- $\epsilon_{ij} \sim N(0, \sigma_\epsilon^2)$ - measurement error
- α - covariate effect
- T_{ik} - covariate indicator

Sialin DataSet (Janssen Pharmaceutica):

- Two groups of mice: wild-type (WT) and knock-out (KO)
- Expression measured on Day 18
- Big changes in gene expression are expected
- Sample size: 6 WT, 6 KO
- ≈ 16000 genes

Number of significant tests

Testing for variance component and treatment effect

		σ_b		
		NS	S	Total
Trt	NS	14910	869	15779
	S	572	44	616
Total		15482	913	16395

Significant Treatment and ICC

Probe Level
Analysis

Khamiakova et
al.

Introduction

Affymetrix
Microarray
Probe Level
Data


Methods

LMM

Dataset

Results and
Discussion

Simulation
Study



Probe Level Analysis

Khamiakova et al.

9 / 17

Significant Treatment, Non-Significant ICC

Probe Level
Analysis

Khiamiakova et
al.

Introduction

Affymetrix
Microarray
Probe Level
Data


Methods

LMM

Dataset

Results and
Discussion

Simulation
Study



Probe Level Analysis

Khiamiakova et al.

10 / 17

Non-Significant Treatment, Significant ICC

Probe Level
Analysis

Khamiakova et
al.

Introduction


Affymetrix
Microarray
Probe Level
Data

Methods

LMM
Dataset

Results and
Discussion

Simulation
Study



Probe Level Analysis

Khamiakova et al.

11 / 17

Non-Significant Treatment, Non-Significant ICC

Probe Level
Analysis

Khiamiakova et
al.

Introduction

Affymetrix
Microarray
Probe Level
Data


Methods

LMM

Dataset

Results and
Discussion

Simulation
Study



Probe Level Analysis

Khiamiakova et al.

12 / 17

Simple Analysis

- on probe level data (estimated treatment effect from LMM)
- t-test on summarized data (FARMS and RMA summarized)
- FDR correction at 0.05

Number of significant genes with treatment effect

	Probe Level	FARMS	RMA
LM, t-test	616	553	371

Results

- more significant genes from probe-level analysis
 - due to higher power (?)
 - due to higher error rate (?)
- Simulation Study
 - investigate power and FDR for probe level analysis
 - varied sample size and probeset size
 - varied treatment effect and intra-class correlation

Results

Probe Level
Analysis

Khmiakova et
al.

Introduction

Affymetrix
Microarray
Probe Level
Data


Methods

LMM

Dataset

Results and
Discussion

Simulation
Study



Analysis of DE

- Probe-level analysis has higher power
- FDR is controlled for all methods

Follow-Up

- 1 Power of other statistical tests on probe level (e.g. SAM)
- 2 Power of the probe level analysis compared to the filtered results
- 3 Extensions of the LMM

Acknowledgments


IWT project

J&J PRD

Hinrich Goehlmann
Willem Talloen
An De Bondt
Luc Bijmens

I-Biostat:

Adetayo Kasim
Ziv Shkedy
Dan Lin
Suzy Van Sanden



Institute of Bioinformatics (Johannes Kepler University):

Sepp Hochreiter
Djork-Arne Clevert

Probe Level
Analysis

Khaniakova et
al.

Introduction

Affymetrix
Microarray
Probe Level
Data

Methods

LMM
Dataset

Results and
Discussion

Simulation
Study